



## Resorbable Haemostatic Sponge

Absorbable Gelatin Sponge USP

# HEMOSPONGE®



ISO 13485:2016  
CERTIFIED COMPANY



INDIAN FDA  
APPROVED



e- Catalogue

# POWER OF 3D MATRIX



## **FAST AND EFFECTIVE HAEMOSTAT**

Stops oozing blood in 2-4 minutes



## **ACTIVATES NATURAL CLOTTING PROCESS ONLY**

Provides surface to auto activate platelets



## **HIGHLY ABSORBENT**

Absorbs and holds the blood upto 45 times the weight of sponge



## **HIGH SURFACE AREA AND UNIFORM POROSITY MATRIX**

Faster clotting and better clot holding capacity



## **100% RESORBABLE**

Completely resorbed in body within four weeks



## **PLIABLE**

Product can be rolled, cut, shaped to fit bleeding site



## **BIO-COMPATIBLE**

Kinder to body tissues (complies ISO 10993 - Biocompatibility Standard)



## **PURE AND SAFE**

Made up of pure pharma grade gelatin only, no chemicals, no drugs



## **pH NEUTRAL**

Can be soaked with thrombin or other drugs for advanced effect



## **ASSURING QUALITY**


Product complies US Pharmacopoeia 42, NF 37 quality standard

# PRODUCT VARIANTS

	<b>Variant</b>	<b>Sponge Size (in mm)</b>	<b>Product reference</b>	<b>Selling Pack sizes</b>
<b>Regular Sizes</b>	STANDARD	80 x 50 x 10	GL851	2 Pcs/Pack, 10 Pcs/Pack 12 Pcs/Pack, 20 Pcs/Pack
	EXPORT	70 x 50 x 10	GL751	2 Pcs/Pack, 10 Pcs/Pack 12 Pcs/Pack, 20 Pcs/Pack
<b>Dental Range</b>	DENTAL SL	10 x 10 x 10	GL111	10 Pcs/Pack, 24 Pcs/Pack, 32 Pcs/Pack, 50 Pcs/Pack
	DENTAL M	14 x 14 x 07	GL1414	10 Pcs/Pack, 50 Pcs/Pack
	DENTAL L	15 x 15 x 15	GL1515	10 Pcs/Pack, 50 Pcs/Pack
	DENTAL XL	20 x 20 x 07	GL202007-D	20 Pcs/Pack, 48 Pcs/Pack
	DENTAL MAX	30 x 30 x 10	GL303010-D	10 Pcs/Pack, 20 Pcs/Pack
<b>ENT Range</b>	ENT SL	20 x 20 x 07	GL202007-E	20 Pcs/Pack, 48 Pcs/Pack
	ENT XL	30 x 30 x 10	GL303010-E	10 Pcs/Pack, 20 Pcs/Pack
	EPISTAXIS	80 x 15 x 08	GL10015	10 Pcs/Pack
	NASAL TAMPON	80 x 10 (Dia)	GL8010	10 Pcs/Pack, 20 Pcs/Pack
	WICK SL	15 x 09	GL1509	10 Pcs/Pack, 20 Pcs/Pack
	WICK XL	24 x 09	GL2409	10 Pcs/Pack, 20 Pcs/Pack
<b>Tampon</b>	TAMPON	80 x 30 (Dia.)	GL8030	2 Pcs/Pack, 6 Pcs/Pack, 10 Pcs/Pack
	TAMPON – SL	60 x 30 (Dia.)	GL6030	2 Pcs/Pack, 6 Pcs/Pack, 10 Pcs/Pack
<b>Powder</b>	POWDER	1 Gram	GL001	6 Pcs/Pack
<b>Film Range</b>	SPECIAL - 1	80 x 50 x 01	GL805001	10 Pcs/Pack
	SPECIAL - 2	70 x 50 x 01	GL705001	10 Pcs/Pack
	FILM – 80	80 x 25 x 01	GL80	12 Pcs/Pack
	FILM – 120	120 x 100 x 01	GL120	12 Pcs/Pack
	FILM – 240	240 x 120 x 01	GL240	12 Pcs/Pack
	FILM – 200	200 x 70 x 0.5	GL200	12 Pcs/Pack

**Custom Size:** HEMOSPONGE is also available in many sizes other than mentioned above, please contact on [info@goodwill-lifesciences.com](mailto:info@goodwill-lifesciences.com) for more details.

# HEMOSPONGE REGULAR SIZES



**For all Surgical  
procedures**

- Rectangle shaped sponge with uniform porosity
- Available in wide variety of packing size
- Can be cut, rolled or mould to desired shape to fit any bleeding site



**Cut, Rolled, Shaped**

**General surgeries | Laparoscopic Surgeries | Neurology  
Cardiology | Orthopaedic | Nephrology | Oncological Surgeries  
Gynaecology | Proctology | Oral, Dental and ENT**



# HEMOSPONGE TAMPONS



**Coloproctological  
Surgeries**




**Gynecological  
Surgeries**

- Cylindrically-shaped sponge with hollow region in centre
- Easily inserted with the help of proctoscope
- Absorbs continuously oozing fluids for longer period
- Enhances patient comfort by reducing post-operative pain and irritation
- Expelled easily with first natural bowl movement



**Post -Haemorrhoidal  
pain relief**



# HEMOSPONGE DENTAL RANGE



**For All  
Dental  
Procedures**



**Small Size  
Big Effect**



**Cost  
Effective**



**Replacement  
of Cotton**

- Small size sponges for use in simple extraction to orthodontic and maxillofacial -oral surgery
- No need to cut - easily fit bleeding sockets
- Better than cotton based products - No chance of re-bleeding
- Can be soaked with thrombin for quick haemostasis
- Can be soaked with antibiotic – for local antiseptic effect
- No need to remove after surgery – provides higher patient comfort
- Individual and sterile piece packing – ready to use and cost - effective
- Made up of pure gelatin - Harmless if swallowed

# HEMOSPONGE ENT RANGE



**Special Sizes to  
fit ENT cavities**



- Specially shaped sponges to fit into small surgical cavities
- Widely accepted and effective haemostatic - packing material
- 100% absorbable - avoids patient discomfort associated with removal of non-absorbable packing material.
- Works as a scaffold to support growth of near-by tissue and fastens the healing process of ENT surgeries.

# HEMOSPONGE FILM RANGE

0.5-  
1.0  
mm


Thickness Range



- Useful for surgeries where large size of sponge is impossible to place i.e. in laminectomy
- Designed for use in thoracic, ocular and neurosurgery.
- Available in various sizes with thickness ranging from 0.5 to 1.00 mm



# HEMOSPONGE POWDER



For Uneven Surface

- Highly effective hemostatic for the sites with uneven surface
- Mixed with saline and kneaded to form putty of desired consistency before application
- The putty has hemostatic effect and is applied directly to bleeding surface

# ABOUT COMPANY

An ISO 14644 Cleanrooms Production Line | An ISO 13485:2016 Certified QMS  
CDSCO Approved Class D Manufacturing Unit | In House Quality Control Laboratories



## e-Catalogues



All variants





Dental Range




ENT Range

# GET IN TOUCH


 [www.goodwill-lifesciences.com](http://www.goodwill-lifesciences.com)

 [info@goodwill-lifesciences.com](mailto:info@goodwill-lifesciences.com)  
[sales@goodwill-lifesciences.com](mailto:sales@goodwill-lifesciences.com)

 +91-87586 18110  
+91-97147 44173

---

### Goodwill Lifesciences Private Limited

 283/284/285, Block no. 392/A, Anjani Textile Park - 4,  
Segwa-Paria Road, At- Paria - 394130, Ta. Olpad, Surat, Gujarat, India  
Mfg Lic no.: MFG/MD/2022/000029



Follow us  
for updates