

## Resorbable Haemostatic Sponge/Tampon



# Absorbable Gelatin Sponge USP HEMOSPONGE<sup>®</sup> *Proctology Range*



ISO 13485:2016  
CERTIFIED COMPANY



INDIAN FDA  
APPROVED



e- Catalogue

# HEMOSPONGE<sup>®</sup>

## Tampon



### Enhances patient Comfort

Reduces the irritation of mucosal membrane of anorectal region



### Cylindrically shaped

Easily placed with help of proctoscope at the end of surgery



### Fast and Effective haemostat

Reduces the continuously oozing post-operative bleeding for prolonged period



### Activates natural clotting process only

Provides surface to auto activate platelets



### Highly Absorbent

Absorbs and holds the blood upto 45 times the weight of sponge



### Pure and Safe

Made up of 100% pure parenteral grade gelatin only, no chemicals, no drugs



### Bio-compatible

Non - toxic, non-allergic and Kinder to body tissues  
(complies ISO 10993 - Biocompatibility Standard)




### pH neutral

Can be soaked with thrombin/antibiotics or other drugs for advanced effect




### Assuring Quality

Product complies US Pharmacopoeia 42, NF 37 quality standard



**High Surface Area and Uniform Porosity Matrix**  
Faster clotting and better clot holding capacity, provides physical supports to overall structure



**100% Resorbable / Easily Expelled**  
Completely resorbed in body within four weeks when left in situ,  
Expelled easily within 1 to 2 days with patient's natural bowl movement when left in anal canal

# Product Variants

Variant	Sponge Size (in mm)	Product reference	Selling Pack sizes
Tampon	80 x 30 (Dia.)	Gl8030	10 Pcs/Pack
Endoscopic Tampon	80 x 10 (Dia.)	Gl8010	10 Pcs/Pack
Standard	80 x 50 x 10	Gl851	02 Pcs/Pack, 10 Pcs/Pack

**Custom Size:** HEMOSPONGE is also available in many sizes other than mentioned above. Please contact on [info@goodwill-lifesciences.com](mailto:info@goodwill-lifesciences.com) for more details.

## e-Catalogues



All variants





Dental Range




ENT Range


## Get in Touch

 [www.goodwill-lifesciences.com](http://www.goodwill-lifesciences.com)

 [info@goodwill-lifesciences.com](mailto:info@goodwill-lifesciences.com)  
[sales@goodwill-lifesciences.com](mailto:sales@goodwill-lifesciences.com)

 +91-87586 18110  
+91-97147 44173

### Goodwill Lifesciences Private Limited

 283/284/285, Block no. 392/A, Anjani Textile Park - 4,  
Segwa-Paria Road, At- Paria - 394130, Ta. Olpad, Surat, Gujarat, India  
Mfg Lic no.: MFG/MD/2022/000029



Follow us  
for updates